

Product datasheet for MR205586

Wnt6 (NM_009526) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Wnt6 (NM_009526) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Wnt6
Synonyms:	AA409270; Wnt-6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR205586 representing NM_009526 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTGCCCGGTGCCCTCCGCCTCGGACTGCTGCTGCTCTTGTGCCCGCGCACGTCGATGGAC
TGTGGTGGCCGTGGGCAGCCCTTGGTCATGGATCCTACCAGCATCTGCAGGAAGGCCAGGCGGTGCC
AGGAAGACAGGCCGAGCTGTGCCAGGCGGAGCCGGAAGTAGTGGCAGAGCTTGCCGAGGCCAAGACTG
GGGTTTCGAGAATGTCAGTTCAGTTCGTTCCGTTCCGACGCTGGAAGTCTCCAGCCACAGCAAGCCCTTG
GGCGCTCCTGCAGCAGGACATCCGAGAGACAGCTTTCGTGTTTGCAATCACCGCAGCTGGTGCCAGCCA
CGCGGTCACTCAAGCCTGTTCCATGGGAGAGCTCCTACAGTGTGGTTGTCAGGCACCCCGGGCGGGCA
CCGCCTAGGCCCTCCGGCCTTCTGGGCACTCCTGGACCTCCAGGACCAACTGGCTCTCCAGATGCTAGCG
CAGCCTGGGAGTGGGAGGCTGCGGAGACGATGTGGACTTCGGGGATGAGAAGTCAAGACTCTTTATGGA
TGCGCAGCACAAGCGGGCCGTGGAGATATCCGTGCATTGGTGAAGTGCACAACAACGAGGCGGGCAGG
CTGGCGGTGCGGAGTACACGCGCACCGAGTGTAAAGTGCATGGGCTTTCGGGTTCTGCGCTCTGCGCA
CCTGCTGGCAGAAGCTGCCTCCGTTCCGCGAGGTGGGCGCACGGCTGCTGGAGCGCTTCCACGGCGCCTC
GCGCGTCATGGCACCAACGACGGCAAAGCTCTGCTGCCTGCGG[AT]CCGCACACTCAAGCCTCCCGGA
CGAGCGGATCTCCTCTACGCAGCCGATTACCCGACTTCTGCGCCCAACCGGCGCACGGGTTCCGCCGG
GCACGCGCGGACGCGCCTGCAACAGCAGTGCCTCCGGACCTCAGCGGCTGCGACCTGTTGTGCTGCGGTG
CGGCACCGCCAGGAGAGCGTACAGCTCAGGAGAAGTCTGTGCCGCTTCCACTGGTGTGCGTGGTG
CAATGCCACCGCTGCCGGTGCGCAAGGAAGTACAGCCTGTGCCTC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTAA



[View online »](#)

Protein Sequence: >MR205586 representing NM_009526
Red=Cloning site Green=Tags(s)

MLPPVPSRLGLLLLLLCPAHVDGLWWAVGSPLVMDPTSI CRKARLAGRQAE LCQAEPEVVAELARGARL
 GVRECQFQFRFRWNCSSHSKAFGRVLQDIRETAFVFAITAAGASHAVTQACSMGELLQCGCQAPRGRA
 PPRPSGLLGTPGPPGPTGSPDASAAWEGGCGDDVDFGDEKSR LFMDAQHKRGRGDIRALVQLHNNEAGR
 LAVRSHTRTECKCHL SGSCALRTCWQKLPFFREVGARLLERFHGASRVMTNDGKALLPAXXRTLKPPG
 RADLLYAADSPDFCAPNRR TGSPGTRGRACNSSAPDL SGCDLLCCGRGHRQESVQLEENCLCRFHWCVV
 QCHRCRVRKELSLCL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9005_e06.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_009526


ORF Size: 1095 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_009526.1 , NM_009526.2 , NM_009526.3 , NP_033552.2
RefSeq Size:	2068 bp
RefSeq ORF:	1095 bp
Locus ID:	22420
UniProt ID:	P22727
Cytogenetics:	1 C4
MW:	40.1 kDa
Gene Summary:	Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. May be a signaling molecule which affects the development of discrete regions of tissues. Is likely to signal over only few cell diameters.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR205586