

## Product datasheet for **MR205587L4V**

### Zdhhc9 (BC042618) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Zdhhc9 (BC042618) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Zdhhc9
Synonyms:	6430508G22; 9530098M12Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC042618
ORF Size:	1092 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205587).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">BC042618</a> , <a href="#">AAH42618</a>
RefSeq Size:	2892 bp
RefSeq ORF:	1094 bp
Locus ID:	208884
Cytogenetics:	X A5
Gene Summary:	The ZDHHC9-GOLGA7 complex is a palmitoyltransferase specific for HRAS and NRAS. [UniProtKB/Swiss-Prot Function]



[View online »](#)