

Product datasheet for **MR205611A1V**

Mouse Fibp (BC049284) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Fibp (BC049284) AAV Particle
Tag:	Myc-DDK
Symbol:	Fibp
Synonyms:	2010004G08Rik; 2010005N21Rik; 3010027N18Rik
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR205611 representing BC049284 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGACCAGCGAACTAGACATTTTCGTGGGGAACACGACCCTTATAGATGAAGACGTGTATCGCCTCTGGC
TGGATGGTTACTCAGTGAACGATGCAGTGGCTCTGCGAGTACGCTCCGGAATCTTGAGCAGACGGGAGC
CACCACAGGAGTGCTGCAGAGCGACACCATGGACCACTACCGCACCTTTCACATGCTTGAGCGTCTGCTG
CACGCGCCGCCGAAGCTGCTGCACCACTCATCTTCCAGATTCCCTCCCGACAGACACTCCTCATCG
AGAGGTACTACACCTTTGATGAGGCCTTTGTTCCGGGAGGTCTTGGCAAGAAGCTGTCCAAGGTACCAA
GAAAGACCTGGATGACATCAGCACCAAAACAGGAATTACTCTCAAGAGCTGCCGGAGGCAGTTTGACAAC
TTTAAACGAGTCTTCAAGGTGGTGAAGAAATGCGGGGCTCCCTGGTGGACAACATCCAGCAGCACTTCC
TCCTCTCTGACCGGCTAGCCAGGGATTACGCAGCCATCGTCTTTTGGCAACAACCGCTTTGAAACAGG
AAAGAAAAAGCTGCAGTACCTGAGCTTTGGGGACTTTGCCTTCTGTGCAGAGCTTATGATCCAGAAGTGG
ACCTTTGGAGCCGTCGGTGAGGCCCCCACTGACCCAGACTCTCAGGTGGACGACATGGATGTGGACTTAG
ATAAGGAGTTTCTCAAGACTTGAAGGAGCTCAAGGTTCTCGTGGCTGACAAGGACCTCTGGACCTGCA
TAAGAGCCTGGTGTGCACTGCCCTCCGGGAAAGCTGGGTGTCTCTGAGATGGAACCAACTTCAAG
AATCTGTCTCGGGGGCTGGTGAACGTGGCTGCCAAGCTGACCCACAATAAGGATGTCAGAGACCTATTTG
TGGACCTCGTGGAGAAGTTTGTGAACCTGCCGCTCTGACCACTGGCCACTGAGTGATGTGCGGCTCTT
CCTCAGCCAGTATTCAGCGTCAGTCCACTCCCTGGATGGCTCCGGCACCAAGGCACTCTGGGACCGCTAC
ATGGGCACCTCCGTGGCTGCCTTCTGCGCCTCTATCATGAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR205611 representing BC049284 Red=Cloning site Green=Tags(s) MTSELDIFVGNTTLIDEDVYRLWLDGYSVNDAVALRVRSGILEQTGATTGVLQSDTMDHYRTFHMLERLL HAPPKLLHQLIFQIPPSRQTLLIERYYTFDEAFVREVLGKKLSKGTKKDLDDISTKTGITLKSCRRQFDN FKRVFKVVEEMRGSVDNIQQHFLLSDRLARDYAAIVFFANNRFETGKKKLQYLSFGDFAFCAELMIQNW TLGAVGEAPDPSQVDDMDVDLDKEFLQDLKELKVLVADKDLLDLHKS LVCTALRGKLGVFSEMETNFK NLSRGLVNVAAKLTHNKDVRDLFVDLVEKFVEPCRS DHWPLSDVRLFLSQYSASVHSLDGRHQALWDRY MGLTRGCLRLYHD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	BC049284
ORF Size:	1092 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	BC049284 , AAH49284
RefSeq Size:	1236 bp
RefSeq ORF:	1094 bp
Locus ID:	58249
Cytogenetics:	19 A
MW:	41.9 kDa