

Product datasheet for **MR205621L3V**

Tarbp2 (NM_009319) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Tarbp2 (NM_009319) Mouse Tagged ORF Clone Lentiviral Particle |
| Symbol: | Tarbp2 |
| Synonyms: | Prbp; TRBP |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_009319 |
| ORF Size: | 1098 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR205621). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_009319.2 |
| RefSeq Size: | 1806 bp |
| RefSeq ORF: | 1098 bp |
| Locus ID: | 21357 |
| UniProt ID: | P97473 |
| Cytogenetics: | 15 57.65 cM |



[View online »](#)

Gene Summary:

Required for formation of the RNA induced silencing complex (RISC). Component of the RISC loading complex (RLC), also known as the micro-RNA (miRNA) loading complex (miRLC), which is composed of DICER1, AGO2 and TARBP2. Within the RLC/miRLC, DICER1 and TARBP2 are required to process precursor miRNAs (pre-miRNAs) to mature miRNAs and then load them onto AGO2. AGO2 bound to the mature miRNA constitutes the minimal RISC and may subsequently dissociate from DICER1 and TARBP2. May also play a role in the production of short interfering RNAs (siRNAs) from double-stranded RNA (dsRNA) by DICER1 (By similarity). Binds in vitro to the PRM1 3' UTR. Seems to act as a repressor of translation (PubMed:8649414).[UniProtKB/Swiss-Prot Function]