

Product datasheet for **MR205622L4V**

Rrs1 (NM_021511) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rrs1 (NM_021511) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rrs1
Synonyms:	5730466A07Rik; AA415037; D1ErtD701e; Rrr
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_021511
ORF Size:	1095 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205622).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_021511.2 , NP_067486.2
RefSeq Size:	2058 bp
RefSeq ORF:	1098 bp
Locus ID:	59014
UniProt ID:	Q9CYH6
Cytogenetics:	1 2.08 cM
Gene Summary:	Involved in ribosomal large subunit assembly. May regulate the localization of the 5S RNP/5S ribonucleoprotein particle to the nucleolus.[UniProtKB/Swiss-Prot Function]


[View online »](#)