

Product datasheet for **MR205641L3V**

Kcnab2 (NM_010598) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Kcnab2 (NM_010598) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Kcnab2
Synonyms:	F5; l2rf5; Kcnb3; kv-beta-2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_010598
ORF Size:	1101 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205641).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_010598.2
RefSeq Size:	3571 bp
RefSeq ORF:	1104 bp
Locus ID:	16498
UniProt ID:	P62482
Cytogenetics:	4 83.08 cM



[View online »](#)

Gene Summary:

Cytoplasmic potassium channel subunit that modulates the characteristics of the channel-forming alpha-subunits (PubMed:8576199). Contributes to the regulation of nerve signaling, and prevents neuronal hyperexcitability (PubMed:11825900, PubMed:21209188). Promotes expression of the pore-forming alpha subunits at the cell membrane, and thereby increases channel activity (By similarity). Promotes potassium channel closure via a mechanism that does not involve physical obstruction of the channel pore (PubMed:8576199). Modulates the functional properties of KCNA4 (By similarity). Modulates the functional properties of KCNA5 (PubMed:8576199). Enhances KCNB2 channel activity (PubMed:8824288). Modulates the functional properties of KCNA5 (PubMed:8576199). Binds NADPH and has NADPH-dependent aldoketoreductase activity (By similarity). Has broad substrate specificity and can catalyze the reduction of methylglyoxal, 9,10-phenanthrenequinone, prostaglandin J2, 4-nitrobenzaldehyde, 4-nitroacetophenone and 4-oxo-trans-2-nonenal (in vitro) (By similarity). [UniProtKB/Swiss-Prot Function]