

## Product datasheet for **MR205652L3V**

### Dusp1 (NM\_013642) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Dusp1 (NM_013642) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dusp1
Synonyms:	3CH134; erp; mkp-1; MKP1; Ptpn16; U19515
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013642
ORF Size:	1104 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205652).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_013642.1</a> , <a href="#">NP_038670.1</a>
RefSeq Size:	1952 bp
RefSeq ORF:	1104 bp
Locus ID:	19252
UniProt ID:	<a href="#">P28563</a>
Cytogenetics:	17 13.28 cM
Gene Summary:	Dual specificity phosphatase that dephosphorylates MAP kinase MAPK1/ERK2 on both 'Thr-183' and 'Tyr-185', regulating its activity during the meiotic cell cycle.[UniProtKB/Swiss-Prot Function]



[View online »](#)