

Product datasheet for **MR205655L4V**

Elavl3 (NM_010487) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Elavl3
Synonyms:	2600009P04Rik; Huc; mHuC; PLE21
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_010487
ORF Size:	1101 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR205655).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_010487.1
RefSeq Size:	5027 bp
RefSeq ORF:	1104 bp
Locus ID:	15571
UniProt ID:	Q60900
Cytogenetics:	9 8.07 cM



Gene Summary:

Binds to AU-rich sequences (AREs) of target mRNAs, including VEGF mRNA. May also bind poly-A tracts via RRM 3. May be involved in neuronal differentiation and maintenance.
[UniProtKB/Swiss-Prot Function]