

Product datasheet for **MR205657L4V**

Rbm39 (BC030493) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rbm39 (BC030493) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rbm39
Synonyms:	1500012C14Rik; 2310040E03Rik; B330012G18Rik; C79248; caper; R75070; Rnpc2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC030493
ORF Size:	1101 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205657).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC030493 , AAH30493
RefSeq Size:	2413 bp
RefSeq ORF:	1103 bp
Locus ID:	170791
Cytogenetics:	2 H1
Gene Summary:	Transcriptional coactivator for steroid nuclear receptors ESR1/ER-alpha and ESR2/ER-beta, and JUN/AP-1. May be involved in pre-mRNA splicing process.[UniProtKB/Swiss-Prot Function]



[View online »](#)