

Product datasheet for **MR205671L3V**

Fbxo28 (NM_175127) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fbxo28 (NM_175127) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fbxo28
Synonyms:	4833428J17Rik; 5730505P19Rik; D1ErtD578e; Fbx28; mKIAA0483
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_175127
ORF Size:	1104 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205671).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_175127.2 , NP_780336.1
RefSeq Size:	4932 bp
RefSeq ORF:	1107 bp
Locus ID:	67948
UniProt ID:	Q8BIG4
Cytogenetics:	1 84.93 cM
Gene Summary:	Probably recognizes and binds to some phosphorylated proteins and promotes their ubiquitination and degradation.[UniProtKB/Swiss-Prot Function]



[View online »](#)