

## Product datasheet for **MR205694**

### Gdap111 (NM\_144891) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Gdap111 (NM\_144891) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Gdap111  
**Synonyms:** Gdap1l  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR205694 representing NM\_144891  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGCGACCCCAACAACCTGACCCCACTGCAGCTGGTGGCCATCTCGGCACTGGAGAGCGATG  
CAGCCAAGCCGGTGGAGACCCCGACGCGCCGAGGCATCCAGCCCGCCATTGGCCAAAGGAGAGCCT  
GGTTCTGTACCACTGGACCCAGTCTTCAGCTCGCAGAAGGTGCGGCTGGTGATCGCTGAGAAGGGCCTG  
GCCTGTGAAGAGAGAGATGTGAGCCTGCCTCAGAGCGAGCATAAGGAACCGTGGTTCATGCGGCTCAACC  
TGGGTGAGGAGGTCCAGTCATTATCCACAGAGACAACATTATCAGTGACTATGACCAGATCATTGACTA  
CGTGGAACGCACCTTCACGGGAGAGCATGTGGTGGCCTTGATGCCCGAGGCGGGCAGCCCAACACGCC  
CGAGTGCTGCAGTACCGGAGCTGCTGGATGCGCTGCCATGGACGCTACACACACGGTGCATCCTGC  
ACCCTGAGCTCACCACAGATTCCATGATCCCAAGTATGCCACGGCTGAGATCCGAGACATTTAGCCAA  
CGCCACCACAGACCTCATGAAACTGGACCACGAGGAGGAACCACAGCTCTCTGAGCCCTACCTTTCTAAA  
CAAAGAAGCTTATGGCCAAGATCCTGGAGCAGATGACGTGAGCTACCTGAAGAAGATCCTCGGGGAGC  
TGGCTATGGTGTGACCAGATTGAGGCGGAGCTGGAGAAGAGGAAACTGGAGAATGAGGGGCAGACGTG  
CGAGCTGTGGCTCTGTGGCTGTGCCTTCACCCTGGCGGATGTCTCTCTGGGAGCCACCCTGCACCCGCTC  
AAGTTCCTGGGACTGTCCAAGAAATACTGGGAAGACGGCAGCCGCCCCAACCTGCAGTCTCTTTGAGA  
GGGTCCAGAGACGCTTTGCCTTCGGAAGTCTGGGCGACATCCACCACTCTGCTGTACGCTGTTAT  
CCCCAACGCGTTCGACTGGTCAAACGGAACCGCCATCCTTCTTTGGGGCATCCTTCTCATGGGCTCC  
TTGGTGGGATGGGCTACTTTGCCTACTGGTACCTCAAGAAAAATACATC

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC  
TGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >MR205694 representing NM\_144891  
Red=Cloning site Green=Tags(s)


MATPNNLTPTNCSWWPISALESDAAKPVETPDAPASSPAHWPKESLVL YHWTQSFSSQKVRLVIAEKGL  
 ACEERDVSLPQSEHKPEWFMRLNLGEEVPVIIHRDNIISDYDQIIDYVERTFTGEHVVALMPEAGSPQHA  
 RVLQYRELLDALPMDAYTHGCILHPELTTDSMIPKYATAEIRRHLANATDLMKLDHEEEPQLSEPYLSK  
 QKKLMAKILEHDDVSYLKKILGELAMVLDQIEAELEKRKLENEGQTCELWLCGCAFTLADVLLGATLHRL  
 KFLGLSKKYWEDGSRPNLQSFFERVQRRAFRKVLGDIHTLLSAVIPNAFRLVKRKPSPFFGASFLMGS  
 LGGMGYFAYWYLKKKYI

SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9019\\_d09.zip](https://cdn.origene.com/chromatograms/mm9019_d09.zip)

**Restriction Sites:** SgfI-RsrII

**Cloning Scheme:**



**ACCN:** NM\_144891

**ORF Size:** 1101 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_144891.2](#), [NP\\_659140.2](#)

**RefSeq Size:** 2580 bp


**RefSeq ORF:** 1104 bp

**Locus ID:** 228858

**Cytogenetics:** 2 H3

**MW:** 42.4 kDa

**Product images:**



Circular map for MR205694