

Product datasheet for **MR205701L3V**

Arpc1a (NM_019767) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Arpc1a (NM_019767) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Arpc1a
Synonyms:	41kDa; 0610010H08Rik; 1110030K07Rik; AA407347; Sid32; Sid329
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_019767
ORF Size:	1113 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205701).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_019767.2 , NP_062741.1
RefSeq Size:	1588 bp
RefSeq ORF:	1113 bp
Locus ID:	56443
UniProt ID:	Q9R0Q6
Cytogenetics:	5 G2



[View online »](#)

Gene Summary:

Probably functions as component of the Arp2/3 complex which is involved in regulation of actin polymerization and together with an activating nucleation-promoting factor (NPF) mediates the formation of branched actin networks (By similarity). In addition to its role in the cytoplasmic cytoskeleton, the Arp2/3 complex also promotes actin polymerization in the nucleus, thereby regulating gene transcription and repair of damaged DNA (By similarity). [UniProtKB/Swiss-Prot Function]