

## Product datasheet for **MR205702L4V**

### Wnt1 (NM\_021279) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Wnt1 (NM_021279) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Wnt1
Synonyms:	Int-1; sw; swaying; Wnt-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_021279
ORF Size:	1110 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205702).
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_021279.4</a>
RefSeq Size:	2371 bp
RefSeq ORF:	1113 bp



[View online >](#)

**Locus ID:** 22408

**UniProt ID:** [P04426](#)

**Cytogenetics:** 15 54.65 cM

**Gene Summary:** Ligand for members of the frizzled family of seven transmembrane receptors. Acts in the canonical Wnt signaling pathway by promoting beta-catenin-dependent transcriptional activation (By similarity). In some developmental processes, is also a ligand for the coreceptor RYK, thus triggering Wnt signaling (PubMed:15454084, PubMed:16116452). Plays an essential role in the development of the embryonic brain and central nervous system (CNS) (PubMed:2202907, PubMed:16116452). Has a role in osteoblast function, bone development and bone homeostasis (By similarity).[UniProtKB/Swiss-Prot Function]