

Product datasheet for **MR205708L3V**

Snrnp40 (NM_025645) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Snrnp40 (NM_025645) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Snrnp40
Synonyms:	0610009C03Rik; Prp8bp; Wdr57
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_025645
ORF Size:	1077 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205708).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025645.2 , NP_079921.2
RefSeq Size:	1613 bp
RefSeq ORF:	1077 bp
Locus ID:	66585
UniProt ID:	Q6PE01
Cytogenetics:	4 D2.2
Gene Summary:	Required for pre-mRNA splicing as component of the activated spliceosome. Component of the U5 small nuclear ribonucleoprotein (snRNP) complex and the U4/U6-U5 tri-snRNP complex, building blocks of the spliceosome.[UniProtKB/Swiss-Prot Function]



[View online »](#)