

Product datasheet for **MR205726L3V**

St13 (NM_133726) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	St13 (NM_133726) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	St13
Synonyms:	1110007I03Rik; 3110002K08Rik; AW555194; HIP; HOP; HSPABP; HSPABP1; p48; PRO0786; SNC6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_133726
ORF Size:	1116 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205726).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_133726.2
RefSeq Size:	1622 bp
RefSeq ORF:	1116 bp
Locus ID:	70356
UniProt ID:	Q99L47
Cytogenetics:	15 E1



[View online »](#)

Gene Summary:

One HIP oligomer binds the ATPase domains of at least two HSC70 molecules dependent on activation of the HSC70 ATPase by HSP40. Stabilizes the ADP state of HSC70 that has a high affinity for substrate protein. Through its own chaperone activity, it may contribute to the interaction of HSC70 with various target proteins (By similarity).[UniProtKB/Swiss-Prot Function]