

## Product datasheet for **MR205736L4V**

### Parva (NM\_020606) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Parva (NM_020606) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Parva
Synonyms:	2010012A22Rik; 5430400F08Rik; Actp; AI225929; AU042898; CH-ILKBP; Parvin
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_020606
ORF Size:	1119 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205736).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_020606.5</a> , <a href="#">NP_065631.3</a>
RefSeq Size:	4450 bp
RefSeq ORF:	1119 bp
Locus ID:	57342
UniProt ID:	<a href="#">Q9EPC1</a>
Cytogenetics:	7 F1



[View online »](#)

**Gene Summary:**

Plays a role in the reorganization of the actin cytoskeleton, formation of lamellipodia and ciliogenesis. Plays a role in the establishment of cell polarity, cell adhesion, cell spreading, and directed cell migration. Plays a role in sarcomere organization and in smooth muscle cell contraction. Required for normal development of the embryonic cardiovascular system, and for normal septation of the heart outflow tract. Plays a role in sprouting angiogenesis and is required for normal adhesion of vascular smooth muscle cells to endothelial cells during blood vessel development.[UniProtKB/Swiss-Prot Function]