

Product datasheet for **MR205740L3V**

Akap10 (BC054105) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Akap10 (BC054105) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Akap10
Synonyms:	1500031L16Rik; B130049N18Rik; D-AKAP-2; D-AKAP2; PRKA10
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC054105
ORF Size:	1116 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205740).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC054105.1
RefSeq Size:	2836 bp
RefSeq ORF:	1118 bp
Locus ID:	56697
Cytogenetics:	11 B2



[View online »](#)

Gene Summary:

This gene encodes a member of A-kinase anchoring proteins (AKAPs), a family of functionally related proteins that target protein kinase A to discrete locations within the cell. The encoded protein is localized to mitochondria and interacts with both the type I and type II regulatory subunits of PKA. It has been reported that this protein is important for maintaining heart rate and myocardial contractility through its targeting of protein kinase A. In mouse, defects of this gene lead to cardiac arrhythmias and premature death. In humans, polymorphisms in this gene may be associated with increased risk of arrhythmias and sudden cardiac death. [provided by RefSeq, May 2013]