

Product datasheet for **MR205769A1V**

Mouse Hnrnpa1 (BC089340) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Hnrnpa1 (BC089340) AAV Particle
Tag:	Myc-DDK
Symbol:	Hnrnpa1
Synonyms:	Hdp; HDP-1; hnrnp-A1; hnRNP A1; Hnrpa1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR205769 representing BC089340 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTCTAAGTCCGAGTCTCCAAGGAGCCAGAACAGCTGCGGAAGCTCTTCATCGGAGGGCTGAGCTTCG
AAACAACCGACGAGAGTCTGAGGAGCCATTTTGAGCAATGGGGAACACTAACAGACTGTGTGGTAATGAG
AGATCCAAACACCAAGAGATCCAGGGGCTTTGGGTTTGTACATATGCCACTGTGGAAGAAGTGGATGCT
GCCATGAATGCAAGACCACACAAGGTGGATGGAAGAGTTGTGGAACCTAAGAGAGCTGTCTCAAGAGAAG
ATTCTCAGCGACCAGGTGCCACTTAACTGTGAAAAAGATCTTTGTTGGTGGTATTAAGAAGACACTGA
AGAACATCACCTACGAGATTATTTTGAGCAGTATGGGAAGATTGAAGTGATAGAAATTATGACTGACAGA
GGCAGTGGGAAAAAGAGGGGCTTTGCTTTTGTACCTTTGATGACCATGACTCTGTGGATAAGATTGTTA
TTCAGAAATACCATACTGTGAATGGCCACAACCTGTGAAGTAAGAAAGGCTCTGTGCAAGCAAGAGATGGC
TAGTGCTTCATCCAGTCAGAGAGGTGCGAGTGGTTCTGGAACTTTGGTGGTGGTCTGGAGGCGGTTTT
GGTGGCAATGACAATTTTGGTCGAGGAGGGAACCTCAGTGGTCTGGTGGCTTTGGTGGCAGCCGTGGT
GTGGTGGATATGGTGGCAGTGGGATGGCTAATGGATTGGCAATGATGGTGGTTATGGAGGAGGCAG
CCCTGGTTACTCTGGAGGAAGCAGAGGCTATGGAAGTGGTGGACAGGTTATGGAACCAAGGCGAGTGGC
TATGGCGGGAGTGGCAGCTATGACAGCTATAACAACGAGGAGGCGGAGGCGGCTTTGGCGGTGGTAGTG
GAAGCAATTTTGGAGGTGGTGAAGCTACAATGATTTTGGCAATTACAACAATCAGTCTTCCAATTTTGG
GCCGATGAAGGGAGGAACTTTGGAGGAGGAGCTCTGGCCCTTATGGTGGTGGAGGCCAGTACTTTGCT
AAACCACGGAACCAAGGTGGCTATGGCGGTTCCAGCAGCAGCAGTAGCTATGGCAGTGGCAGGAGTTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR205769 representing BC089340 Red=Cloning site Green=Tags(s) MSKSESPKEPEQLRKLFIGGLSFETTTDESLRSHFEQWGTLTDCVVMRDPNPKRSRGFGFVTYATVEEVDA AMNARPHKVDGRVVEPKRAVSREDSQRPGAHLTVKKIFVGGIKEDTEEHHLRDYFEQYGKIEVIEIMTDR GSGKKRGFAFVTDDHDSVDKIVIQYHTVNGHNCEVRKALSKQEMASASSSQGRSGSGNFGGGRGGGF GGNDNFGRRGNFSGRGGFGGSRGGGGYGGSGDGYNGFGNDGGYGGGSPGYSGGSRGYSGGGQYGNQGS YGGSGSYDSYNNGGGGGGFGGGSGSNFGGGGSYNDFGNYNQSSNFGPMKGGNFGRSSGPYGGGGQYFA KPRNQGGYGGSSSSSYGSGRRF TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	BC089340
ORF Size:	1119 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	BC089340.1
RefSeq Size:	1891 bp
RefSeq ORF:	1121 bp
Locus ID:	15382
Cytogenetics:	15 58.58 cM
MW:	38.8 kDa