

## Product datasheet for MR205769L3

### Hnrpa1 (BC089340) Mouse Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Hnrpa1 (BC089340) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Hnrpa1
Synonyms:	Hdp; HDP-1; hnrnp-A1; hnRNP A1; Hnrpa1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205769).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.

ACCN:	BC089340
ORF Size:	1119 bp



[View online »](#)

