

Product datasheet for **MR205769L3V**

Hnrpa1 (BC089340) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Hnrpa1 (BC089340) Mouse Tagged ORF Clone Lentiviral Particle |
| Symbol: | Hnrpa1 |
| Synonyms: | Hdp; HDP-1; hnrnp-A1; hnRNP A1; Hnrpa1 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | BC089340 |
| ORF Size: | 1119 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR205769). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | BC089340.1 |
| RefSeq Size: | 1891 bp |
| RefSeq ORF: | 1121 bp |
| Locus ID: | 15382 |
| Cytogenetics: | 15 58.58 cM |
| Gene Summary: | Involved in the packaging of pre-mRNA into hnRNP particles, transport of poly(A) mRNA from the nucleus to the cytoplasm and may modulate splice site selection (By similarity). May bind to specific miRNA hairpins (By similarity).[UniProtKB/Swiss-Prot Function] |



[View online »](#)