

Product datasheet for **MR205769L4V**

Hnrpa1 (BC089340) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hnrpa1 (BC089340) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hnrpa1
Synonyms:	Hdp; HDP-1; hnrnp-A1; hnRNP A1; Hnrpa1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC089340
ORF Size:	1119 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205769).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC089340.1
RefSeq Size:	1891 bp
RefSeq ORF:	1121 bp
Locus ID:	15382
Cytogenetics:	15 58.58 cM
Gene Summary:	Involved in the packaging of pre-mRNA into hnRNP particles, transport of poly(A) mRNA from the nucleus to the cytoplasm and may modulate splice site selection (By similarity). May bind to specific miRNA hairpins (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)