

Product datasheet for MR205777

Nfe2 (NM_008685) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nfe2 (NM_008685) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Nfe2
Synonyms:	NF-E2; NF-E2/P45; p45; p45nf-e2; p45NFE2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR205777 representing NM_008685 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCCCGTGTCTCTCAGCAGAACAGGAACAGGTTATCACAGCTGCCTGTTGGGGAGCTTGGAGAGA
TGGAACTGACTTGGCAAGAGATCATGTCCACTGAGCTGCAGGGTCTAAATGTTCCAAGTGAGACATC
TTTTGAGCCTCAAGCACCCACCCATACCTGGGCCACTGCCACCTCCAACATATTGCCCTGTTCATT
CATCCAGATGCAGGCTTCTCCCTCCCCACCATCTTATGAGCTCCAGCATCTACTCCCATGTCCAG
AACTACCATACTCTATGGTAATGTAGCCATACCAGTGTCAAAGCCACTTACCCTTCAGGCCTGCTCAA
TGAGCCCCTCCCAGACCACTTAGCTCTCCTGGACATTGGGCTGCCAGTGGGGCAACCAAGCCCCAAGAA
GACCCAGAATCTGACTCAGGATTATCCCTCAACTACAGTGATGCAGAATCTCTTGAGCTAGAGGGTATGG
AGGCTGGCAGGCGGAGGAGCGAGTACGCGGACATGTACCCAGTGGAGTATCCTTACTCACTTATGCCAA
TTCTTTGGCCATCCCAACTATACTCTTCCACCCACTGAGACACCTTGGCCTTAGAGTCATCTCCGGT
CCAGTTCGGGCTAAGCCTGCTGTCCGTGGGGAGGCAGGGAGTCCGGACGAGCGGCGAGCCCTGGCCATGA
AGATTCCTTTCCCTACGGACAAGATAGTTAACTTGCCGGTAGATGACTTTAATGAGTTGTTGGCACAGTA
TCCGCTAACGGAGAGCCAGCTGGCTCTAGTTCGGGACATCCGTCGACGGGCAAGAACAAGTGGCAGCC
CAAACTGTCGCAAGAGAAAACCTGAAAACATTGTGCAGCTGGAGCGAGAGCTGGAGCGGCTGAGCAGTG
AAAGGGAGCGGCTTCTCAGAGCCCGAGGGGAGGCTGACCCCACTCTGGAGGTCATGCGCCAACAGCTGGC
AGAGCTGTACCATGATATTTTCCAGCATCTTCGGGATGAATCTGGCAACAGTTACTACCCAGAGGAATAT
GTAAGTGAACAGGCTGCTGATGGTCCATCTTTCTGGTACCCCGTGAACCAAAATGGAGGCTACAGAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTAA



[View online »](#)

Protein Sequence: >MR205777 representing NM_008685
 Red=Cloning site Green=Tags(s)

MPPCPPQQRNRLSQLPVGELGEMELTWQEIMSITELQGLNVPSETSFEPQAPTPYGPLPPPTYCPCSI
 HPDAGFSLPPSYELPASTPHVPELPYSYGNVAIPVSKPLTSLGLLNEPLPDHLALLDIGLPVGQPKPQE
 DPESDSSLNYSDAESLELEGMAGRRRSEYADMYVPEYPSLMPNSLAHPNYLPPTEPLALESSSG
 PVRAKPAVRGEAGSRDERRALAMKIPFPDKIVNLPVDDFNELLAQYPLTESQLALVRDIRRRGKKNKVA
 QNCRKRKLETIVQLERELERLSSERLLRARGEADRTLEVMRQQLAELYHDFQHLRDESGNSYSPEEY
 VLQQAADGAIFLVPRGTMKMEATD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9039_a06.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



* The last codon before the Stop codon of the ORF

ACCN: NM_008685

ORF Size: 1119 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_008685.3](#), [NP_032711.2](#)

RefSeq Size: 1793 bp

RefSeq ORF: 1122 bp

Locus ID: 18022


UniProt ID: [Q07279](#)

Cytogenetics: 15 58.62 cM

MW: 42 kDa

Gene Summary: Component of the NF-E2 complex essential for regulating erythroid and megakaryocytic maturation and differentiation. Binds to the hypersensitive site 2 (HS2) of the beta-globin control region (LCR). This subunit (NFE2) recognizes the TCAT/C sequence of the AP-1-like core palindrome present in a number of erythroid and megakaryocytic gene promoters. Requires MAFK or other small MAF proteins for binding to the NF-E2 motif. May play a role in all aspects of hemoglobin production from globin and heme synthesis to procurement of iron. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR205777