

## Product datasheet for **MR205790L4V**

### Cxcr5 (NM\_007551) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Cxcr5 (NM_007551) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cxcr5
Synonyms:	Blr1; CXC-R5; CXCR-5; Gpcr6; MDR15
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_007551
ORF Size:	1122 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205790).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_007551.2</a> , <a href="#">NP_031577.2</a>
RefSeq Size:	2636 bp
RefSeq ORF:	1125 bp
Locus ID:	12145
UniProt ID:	<a href="#">Q04683</a>
Cytogenetics:	9 24.84 cM
Gene Summary:	Cytokine receptor that binds to B-lymphocyte chemoattractant (BLC). Involved in B-cell migration into B-cell follicles of spleen and Peyer patches but not into those of mesenteric or peripheral lymph nodes.[UniProtKB/Swiss-Prot Function]



[View online »](#)