

Product datasheet for **MR205812L3V**

Fbxl8 (NM_015821) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Fbxl8
Synonyms:	FBL8
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_015821
ORF Size:	1122 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR205812).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_015821.1](#), [NP_056636.1](#)

RefSeq Size: 1901 bp

RefSeq ORF: 1125 bp

Locus ID: 50788

UniProt ID: [Q8CIG9](#)

Cytogenetics: 8 D3



Gene Summary:

Substrate-recognition component of the SCF (SKP1-CUL1-F-box protein)-type E3 ubiquitin ligase complex.[UniProtKB/Swiss-Prot Function]