

Product datasheet for **MR205822L3V**

Serpinb5 (NM_009257) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Serpinb5 (NM_009257) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Serpinb5
Synonyms:	1110036M19Rik; A1462524; A1646751; Maspin; ovalbumin; PI-5; Spi7
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_009257
ORF Size:	1128 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205822).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_009257.3
RefSeq Size:	2565 bp
RefSeq ORF:	1128 bp
Locus ID:	20724
UniProt ID:	P70124
Cytogenetics:	1 E2.1



[View online »](#)

Gene Summary:

Tumor suppressor. It blocks the growth, invasion, and metastatic properties of mammary tumors. As it does not undergo the S (stressed) to R (relaxed) conformational transition characteristic of active serpins, it exhibits no serine protease inhibitory activity (By similarity). [UniProtKB/Swiss-Prot Function]