

Product datasheet for **MR205837L3V**

Tpst2 (NM_009419) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tpst2 (NM_009419) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tpst2
Synonyms:	AI448750; D5Ucla3; grm; grt; Tango13b
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_009419
ORF Size:	1131 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205837).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_009419.2 , NP_033445.1
RefSeq Size:	1830 bp
RefSeq ORF:	1173 bp
Locus ID:	22022
UniProt ID:	O88856
Cytogenetics:	5 54.65 cM
Gene Summary:	Catalyzes the O-sulfation of tyrosine residues within acidic motifs of polypeptides, using 3'-phosphoadenylyl sulfate (PAPS) as cosubstrate.[UniProtKB/Swiss-Prot Function]



[View online »](#)