

Product datasheet for **MR205838A1V**

Mouse Birc3 (BC011338) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Birc3 (BC011338) AAV Particle
Tag:	Myc-DDK
Symbol:	Birc3
Synonyms:	IAP1, MIAP1, Birc2, cIAP1, cIAP-1, MIHB, HIAP2, Api2, IAP2, MIHC, MIAP2, RNF49, cIAP2, cIAP-2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR205838 representing BC011338 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAACATGGTTCAAGACAGTGCCTTTCTAGCCAAGCTGATGAAGAGTGCTGACACCTTTGAGTTGAAGT
ATGACTTTTCTGTGAGCTGTACCGATTGTCCACATATTCAGCTTTTCCAGGGGAGTTCTGTGTGTCAGA
AAGGAGTCTGGCTCGTGTGGCTTTTACTACACTGGTGCAATGACAAGGTCAAGTGCTTCTGTGTGGC
CTAATGCTAGACAACCTGGAAACAAGGGGCGAGTCCCATGGAGAAGCACAGAAAGTTGTACCCAGCTGCA
ACTTTGTACAGACTTTGAATCCAGCCAACAGTCTGGAAGCTAGTCCTCGGCCTTCTCTTCTTCCACGGC
GATGAGCACCATGCCTTTGAGCTTTGCAAGTTCTGAGAACACTGGCTATTTTCAGTGGCTCTTACTCGAGC
TTTCCCTCAGACCCTGTGAACCTCCGAGCAAAATCAAGATTGTCCTGCTTTGAGCACAAGTCCCTACCACT
TTGCAATGAACACAGAGAAGGCCAGATTACTCACCTATGAAACATGGCCATTGTCTTTTCTGTACCAGC
AAAGCTGGCCAAAGCAGGCTTCTACTACATAGGACCTGGAGATAGAGTGGCCTGCTTTGCGTGCGATGGG
AAACTGAGCAACTGGGAACGTAAGGATGATGCTATGTCAGAGCACCAGAGGCATTTCCAGCTGCCCGT
TCTTAAAGACTTGGGTCAGTCTGCTTCGAGATACACTGTCTCTAACCTGAGCATGCAGACACACGCAGC
CCGTATTAGAACATTCTAAGTGGCCTTCTAGTGCAGTATTCATTCCAGGAAGTTGCAAGTGGGGC
TTTTATTATACAGGACACAGTGATGATGTCAAGTGTGTTTGTGTGATGGTGGGCTGAGGTGCTGGGAAT
CTGGAGATGACCCTGGGTGGAACATGCCAAGTGGTTTCCAAGGTGTGAGTACTTGCTCAGAATCAAAGG
CCAAGAATTTGTAGCCAAGTTCAAGCTGGCTATCCTCATCTACTTGAGCAGCTATTATCTACGTACAGC
TCCCAGAAGATGAGAATGCAGACGCAGCAAGTATGTATAATAATCATAATTCCTGCATTACACTGCTCC
GTTTTGCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



View online »

Protein Sequence:	>MR205838 representing BC011338 Red=Cloning site Green=Tags(s) MNMVQDSAFLAKLMKSADTFELKYDFSCELYRLSTYSAFPGRGVPVSERSLARAGFYTGVDKVKCFCCG LMLDNWKQGGSPMEKHRKLYPSCNFVQTLNPANSLEASPRPSLPSTAMSTMPLSFASSENTGYFSGSYSS FPSDPVNFRANQDCPALSTSPYHFAMNTEKARLLTYETWPLSFLSPAKLAKAGFYIIGPDRVACFACDG KLSNWERKDDAMSEHQHFPSCPFLKDLGQSASRYTVSNLSMQTHAARIRTFSNWPSSALVHSQELASAG FYITGHSDDVKCFCCDGGRLRCWESGDDPWVEHAKWFPCEYLLRIKGQEFVSQVQAGYPHLLLEQLLSTSD SPEDENADAASMYNNHNSCITLLRFA TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	BC011338
ORF Size:	1128 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	BC011338.1
RefSeq Size:	2979 bp
RefSeq ORF:	1130 bp
Locus ID:	11796
Cytogenetics:	9 A1
MW:	109.2 kDa