

Product datasheet for **MR205840A1V**

Mouse Ccdc9b (NM_001001982) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Ccdc9b (NM_001001982) AAV Particle
Tag:	Myc-DDK
Symbol:	Ccdc9b
Synonyms:	Ccdc9b; Gm353
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR205840 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGAACAAAGCCAAGGCCAAAGGATAGTAAGTGAGAAACCCACGACCCAAGCGCGGAGCCCAAGGG
CAAAGGGGACATCTGGAGGAGGCCGGAGCTGGCGTGCACCTTGCAGATGACTTCCAGCTCAGATTCAGC
AGGGAAGGGTGTCTGGAGCCCCGGAGTCCAGGCGTGGTCTCAAGTCCACCTCCCAACAGTCTTTGGGG
CTGTCCCCTGAGGCCAGCTGGGACTACACGCACTGGAACAAGAGAGAGAGCAGATTGACCTGGCCCCGCC
TTGCCCGGCACAGAAACCCACAGGGTGACTGGAGCCCGCTGGGACCTGGACAAAGCCAAGTCCATGCC
ACAGAAGTGCACAAGCCTAGGGACAGGGACCTGGCCATGGGCAGCAGGAAGGGTCTTAGGACCCACCAG
AAGCCTCAGCACCCACCATCATCCCTGGATGGAATGTGCGGGTGGACAATCTGGAAGACCCTCGATGA
CATCCACCATAGGCAGCAACGCCAGGGGAAGGAGAGGCTAAGTGGCAGGGCCCGAAGGTATGACATGAA
AGAAGATGAGATGCTGCAGAGCCAGGAGGGAAGCCAGAGCAGCACGAAGACTCCAGTGCAGAGGAAGAA
TTCGTACAGAAGCGGAATGAAAGGGAGCCAGACAGACAAGAGACCGTGCCTGCAACCAGCCCAATTCCAG
CCTCCCAGAGGGGCTAAAAGGAGAGTCGGGGACCTCTACAGCCAGTTTGGCCCCAGACTCTCCACAACA
TTCTGACTTGGCTCCCCTCGACCTCTCTGTGGGAGGAGCTAGCAGCCGCAAGCCTGGAAGAAAGTACATGT
ATACTCGGTCCAAAGGCTGGGGCTCAGCAAAGTCCAGTGTGATGCCAAATGGTTCTGAGCAGCAGACCC
TGGGGTCACTGACCACAGCCTGGGCTGGAAGTCCACAGTTGCACTGAGTACAGAGGAAAGCAGACAGC
CCCAGAGCCCATTGGAGGACAGGGCTGGCAAGGCTGGGGCCCAAGAGAATCTGACACCAAGAAGCAGGCCC
CCCAGAGGAATTGGCCCAAGGGCAAGAGGCACAAGTAAAGAAGCAGGGCAGGAGGTCTGGCCCTGCAG
GAAGATGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR205840 protein sequence Red=Cloning site Green=Tags(s) MENKAKAKRIVSEKPTTQARSPRAKGTSGGGRSWRAPLQMTSSSDSAGKGVLEPRSPGVVSSPPPPQQSLG LSPEASWDYTQWKQEREQIDLARLARHRNPQGDWSRPWDLKAKSMPQNCHKPRDRDLAMGSRKGPRTHQ KPQHPPSSLDGKCRGGQSGRPSMTSTIGSNAQGKERLTGRARRYDMKEDEMLQSQEGSQSSTKTPSAEEE FVQKRNEREPDRQETVPATSPIPASPEGLKGESGTSTASLAPDSPQHSDLAPLDLSVGGASSRKPGKSTC ILGPKAGAQQSPVSWPNGSEQQLGCTDHQPGLLEVHSCTESQRKAQTPEPMEDRAGKAGAQENLTPRS PRGIGPRARGTSKRSRAGGPGPAGRC TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_001001982
ORF Size:	1128 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001001982.1</u> , <u>NP_001001982.1</u>
RefSeq Size:	4716 bp
RefSeq ORF:	1638 bp
Locus ID:	214239
UniProt ID:	<u>A3KGF9</u>
Cytogenetics:	2 E5
MW:	40.4 kDa