

Product datasheet for MR205847L3

Nif311 (NM_022988) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nif311 (NM_022988) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Nif311
Synonyms:	1110030G24Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205847).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_022988
ORF Size:	1131 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_022988.2
RefSeq Size:	2161 bp
RefSeq ORF:	1131 bp
Locus ID:	65102
UniProt ID:	Q9EQ80
Cytogenetics:	1 C1.3
Gene Summary:	This gene is a member of the NGG1-interacting factor 3-like superfamily of transcriptional regulators and is ubiquitously expressed throughout embryonic development. The encoded protein interacts with a component of the constitutive photomorphogenesis 9 signalosome, and functions as a transcriptional corepressor of genes involved in neuronal differentiation. This gene is highly conserved from bacteria to human. [provided by RefSeq, Sep 2016]

Product images:



Circular map for MR205847L3