

## Product datasheet for **MR205892L3V**

### Ube2q2 (NM\_180600) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ube2q2 (NM_180600) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ube2q2
Synonyms:	3010021M21Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_180600
ORF Size:	1137 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205892).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_180600.1</a> , <a href="#">NP_850931.1</a>
RefSeq Size:	2952 bp
RefSeq ORF:	1137 bp
Locus ID:	109161
UniProt ID:	<a href="#">Q8K2Z8</a>
Cytogenetics:	9 B
Gene Summary:	Accepts ubiquitin from the E1 complex and catalyzes its covalent attachment to other proteins. In vitro catalyzes 'Lys-48'-linked polyubiquitination.[UniProtKB/Swiss-Prot Function]



[View online »](#)