

Product datasheet for **MR205904L4V**

Cdc37 (NM_016742) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cdc37 (NM_016742) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cdc37
Synonyms:	p50; p50Cdc37
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_016742
ORF Size:	1140 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205904).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016742.2
RefSeq Size:	2456 bp
RefSeq ORF:	1140 bp
Locus ID:	12539
UniProt ID:	Q61081
Cytogenetics:	9 A3
Gene Summary:	Co-chaperone that binds to numerous kinases and promotes their interaction with the Hsp90 complex, resulting in stabilization and promotion of their activity. Inhibits HSP90AA1 ATPase activity (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)