

Product datasheet for **MR205909L3V**

Rrh (BC046288) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rrh (BC046288) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rrh
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC046288
ORF Size:	1137 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205909).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC046288 , AAH46288
RefSeq Size:	1533 bp
RefSeq ORF:	1139 bp
Locus ID:	20132
Cytogenetics:	3 59.09 cM
Gene Summary:	May play a role in rpe physiology either by detecting light directly or by monitoring the concentration of retinoids or other photoreceptor-derived compounds.[UniProtKB/Swiss-Prot Function]


[View online »](#)