

## Product datasheet for **MR205909L4V**

### Rrh (BC046288) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | Lentiviral Particles                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Product Name:             | Rrh (BC046288) Mouse Tagged ORF Clone Lentiviral Particle                                                                                                                                                                                                                                                                                                                                                                                          |
| Symbol:                   | Rrh                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Mammalian Cell Selection: | Puromycin                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Vector:                   | pLenti-C-mGFP-P2A-Puro (PS100093)                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Tag:                      | mGFP                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| ACCN:                     | BC046288                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| ORF Size:                 | 1137 bp                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(MR205909).                                                                                                                                                                                                                                                                                                                                                                                     |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.                                                                                                                                                                                                                                                                                                               |
| RefSeq:                   | <a href="#">BC046288</a> , <a href="#">AAH46288</a>                                                                                                                                                                                                                                                                                                                                                                                                |
| RefSeq Size:              | 1533 bp                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| RefSeq ORF:               | 1139 bp                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Locus ID:                 | 20132                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Cytogenetics:             | 3 59.09 cM                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Gene Summary:             | May play a role in rpe physiology either by detecting light directly or by monitoring the concentration of retinoids or other photoreceptor-derived compounds.[UniProtKB/Swiss-Prot Function]                                                                                                                                                                                                                                                      |


[View online »](#)