

Product datasheet for **MR205918L4V**

Hnrnpa3 (NM_198090) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hnrnpa3 (NM_198090) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Hnrnpa3
Synonyms:	2410013L13Rik; 2610209F03Rik; 2610510D13Rik; Hnrpa3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_198090
ORF Size:	1137 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205918).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_198090.2 , NP_932758.1
RefSeq Size:	2639 bp
RefSeq ORF:	1140 bp
Locus ID:	229279
UniProt ID:	Q8BG05
Cytogenetics:	2 C3
Gene Summary:	Plays a role in cytoplasmic trafficking of RNA. Binds to the cis-acting response element, A2RE. May be involved in pre-mRNA splicing (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)