

## Product datasheet for **MR205924L3V**

### Flot2 (NM\_008028) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Flot2 (NM_008028) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Flot2
Synonyms:	AI573412; Esa; reggie-2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_008028
ORF Size:	1140 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205924).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_008028.2</a> , <a href="#">NP_032054.1</a>
RefSeq Size:	2718 bp
RefSeq ORF:	1140 bp
Locus ID:	14252
UniProt ID:	<a href="#">Q60634</a>
Cytogenetics:	11 46.74 cM
Gene Summary:	May act as a scaffolding protein within caveolar membranes, functionally participating in formation of caveolae or caveolae-like vesicles. May be involved in epidermal cell adhesion and epidermal structure and function.[UniProtKB/Swiss-Prot Function]



[View online »](#)