

Product datasheet for **MR205931L3V**

Cers2 (NM_029789) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cers2 (NM_029789) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cers2
Synonyms:	0610013I17Rik; AI225939; Lass2; TRH3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_029789
ORF Size:	1143 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205931).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_029789.1
RefSeq Size:	2048 bp
RefSeq ORF:	1143 bp
Locus ID:	76893
UniProt ID:	Q924Z4
Cytogenetics:	3 F2.1
Gene Summary:	Suppresses the growth of cancer cells. May be involved in sphingolipid synthesis (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)