

Product datasheet for MR205943L3

Rnf13 (NM_001113413) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rnf13 (NM_001113413) Mouse Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Rnf13
Synonyms:	2010001H16Rik; Rz; Rzf
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205943).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_001113413
ORF Size:	1146 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001113413.1
RefSeq Size:	2676 bp
RefSeq ORF:	1146 bp
Locus ID:	24017
UniProt ID:	O54965
Cytogenetics:	3 D
Gene Summary:	This gene encodes a member of the PA-TM-RING family of proteins that contain a protease associated (PA) domain and a RING finger domain separated by a transmembrane (TM) domain. The encoded protein is an E3 ubiquitin ligase localized to the endosomal-lysosomal vesicles and inner nuclear membrane. Mice lacking the encoded protein have impaired learning abilities associated with a decreased synaptic vesicle density and dysregulated SNARE complex assembly. Alternative splicing of this gene results in multiple transcript variants. A pseudogene for this gene has been identified on the X chromosome. [provided by RefSeq, Jan 2015]

Product images:



Circular map for MR205943L3