

Product datasheet for **MR205943L3V**

Rnf13 (NM_001113413) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rnf13 (NM_001113413) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rnf13
Synonyms:	2010001H16Rik; Rz; Rzf
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001113413
ORF Size:	1146 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205943).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001113413.1
RefSeq Size:	2676 bp
RefSeq ORF:	1146 bp
Locus ID:	24017
UniProt ID:	O54965
Cytogenetics:	3 D



[View online »](#)

Gene Summary:

This gene encodes a member of the PA-TM-RING family of proteins that contain a protease associated (PA) domain and a RING finger domain separated by a transmembrane (TM) domain. The encoded protein is an E3 ubiquitin ligase localized to the endosomal-lysosomal vesicles and inner nuclear membrane. Mice lacking the encoded protein have impaired learning abilities associated with a decreased synaptic vesicle density and dysregulated SNARE complex assembly. Alternative splicing of this gene results in multiple transcript variants. A pseudogene for this gene has been identified on the X chromosome. [provided by RefSeq, Jan 2015]