

## Product datasheet for **MR205950L2V**

### **NdrG1 (NM\_008681) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	NdrG1 (NM_008681) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	NdrG1
Synonyms:	CAP43; CMT4D; DRG1; HMSNL; Ndr1; NdrI; NMSL; PROXY1; RTP; TDD5
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_008681
ORF Size:	1185 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205950).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_008681.2</a> , <a href="#">NP_032707.2</a>
RefSeq Size:	2909 bp
RefSeq ORF:	1185 bp
Locus ID:	17988
UniProt ID:	<a href="#">Q62433</a>
Cytogenetics:	15 D2



[View online »](#)

**Gene Summary:**

Stress-responsive protein involved in hormone responses, cell growth, and differentiation. Acts as a tumor suppressor in many cell types. Necessary but not sufficient for p53/TP53-mediated caspase activation and apoptosis. Required for vesicular recycling of CDH1 and TF. May also function in lipid trafficking. Protects cells from spindle disruption damage. Functions in p53/TP53-dependent mitotic spindle checkpoint. Regulates microtubule dynamics and maintains euploidy (By similarity). Has a role in cell trafficking notably of the Schwann cell and is necessary for the maintenance and development of the peripheral nerve myelin sheath. [UniProtKB/Swiss-Prot Function]