

Product datasheet for **MR205964L1V**

Gja1 (NM_010288) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gja1 (NM_010288) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gja1
Synonyms:	AU042049; AW546267; Cnx43; connexin43; Cx43; Cx43alpha1; Gja-1; Npm1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_010288
ORF Size:	1146 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205964).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_010288.2
RefSeq Size:	3071 bp
RefSeq ORF:	1149 bp
Locus ID:	14609
UniProt ID:	P23242
Cytogenetics:	10 28.64 cM



[View online »](#)

Gene Summary:

Gap junction protein that acts as a regulator of bladder capacity. A gap junction consists of a cluster of closely packed pairs of transmembrane channels, the connexons, through which materials of low MW diffuse from one cell to a neighboring cell. Negative regulator of bladder functional capacity: acts by enhancing intercellular electrical and chemical transmission, thus sensitizing bladder muscles to cholinergic neural stimuli and causing them to contract. May play a role in cell growth inhibition through the regulation of NOV expression and localization (PubMed:15181016). Plays an essential role in gap junction communication in the ventricles (PubMed:26403541).[UniProtKB/Swiss-Prot Function]