

Product datasheet for **MR205973L4V**

B4galt6 (NM_019737) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	B4galt6 (NM_019737) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	B4galt6
Synonyms:	AA536803; AU022389
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_019737
ORF Size:	1149 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205973).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_019737.2 , NP_062711.1
RefSeq Size:	5808 bp
RefSeq ORF:	1149 bp
Locus ID:	56386
UniProt ID:	Q9WVK5
Cytogenetics:	18 A2



[View online »](#)

Gene Summary:

Catalyzes the synthesis of lactosylceramide (LacCer) via the transfer of galactose from UDP-galactose to glucosylceramide (GlcCer) (PubMed:23882130, PubMed:30114188). LacCer is the starting point in the biosynthesis of all gangliosides (membrane-bound glycosphingolipids) which play pivotal roles in the CNS including neuronal maturation and axonal and myelin formation (PubMed:30114188).[UniProtKB/Swiss-Prot Function]