

Product datasheet for **MR205990L4V**

Lef1 (BC057543) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lef1 (BC057543) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lef1
Synonyms:	3000002B05; AI451430; Lef-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC057543
ORF Size:	1152 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR205990).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC057543 , AAH57543
RefSeq Size:	2999 bp
RefSeq ORF:	1154 bp
Locus ID:	16842
Cytogenetics:	3 60.78 cM



[View online »](#)

Gene Summary:

Participates in the Wnt signaling pathway. Activates transcription of target genes in the presence of CTNNB1 and EP300. May play a role in hair cell differentiation and follicle morphogenesis. TLE1, TLE2, TLE3 and TLE4 repress transactivation mediated by LEF1 and CTNNB1 (By similarity). Regulates T-cell receptor alpha enhancer function. Binds DNA in a sequence-specific manner. PIASG antagonizes both Wnt-dependent and Wnt-independent activation by LEF1.[UniProtKB/Swiss-Prot Function]