

Product datasheet for **MR206012L3V**

**Lenti ORF particles, MGC:76795 (Myc-DDK-tagged) - Mouse cDNA clone MGC:76795
IMAGE:30111586, 200ul, >10⁷ TU/mL (Gm13152) (BC066875) Mouse Tagged ORF Clone
Lentiviral Particle**

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lenti ORF particles, MGC:76795 (Myc-DDK-tagged) - Mouse cDNA clone MGC:76795 IMAGE:30111586, 200ul, >10 ⁷ TU/mL (Gm13152) (BC066875) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gm13152
Synonyms:	OTTMUSG00000010671
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC066875
ORF Size:	1155 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206012).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC066875 , AAH66875
RefSeq Size:	2891 bp
RefSeq ORF:	1157 bp
Locus ID:	195531
Cytogenetics:	4



[View online »](#)