

Product datasheet for **MR206038A1V**

Mouse Serbp1 (BC013665) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Serbp1 (BC013665) AAV Particle
Tag:	Myc-DDK
Symbol:	Serbp1
Synonyms:	1200009K13Rik; 9330147J08Rik; AL022786; Pairbp1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR206038 representing BC013665 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCCTGGGCACCTACAGGAAGGCTTCGGCTGCGTGGTCACCAACCGATTCGACCAGCTATTTGACGACG
AATCGGACCCTTTGAGGTACTGAAGGCAGCAGAGACAAGAAAAAGAGCCGGCGGGGCGGCGTTGG
GGGCCCCGGGGCCAAGAGCGCGGCTCAGGCCGCGGCCAGACCAACTCCAACGCGGCGGGCAAACAGTTG
CGTAAAGAGTCCCAGAAAGACCGCAAGAACCCTGCCCCCAGCGTCGGCGTGGCCGACAAAAAGGAGG
AGACGCAGCCGCGGTGGCGCTTAAGAAAGAAGGAATAAGGCGAGTTGGAAGAAGACCCGATCAACAAC
ACAGGGTGATGGGAACTAATTGATAGGAGAGCAGAAAGGCGACCAACCCCGTGAAAGAAGATTTGAAAA
CCAATTGAAGAAAAAGGTGAAGGAGGTGAATTTTCAGTTGATAGACCGATTATCGAACGGCCTATCCGAG
GCCGAGGTGGTCTTGAAGGGGTGAGGAGGCCGTGGACGTGGAATGGGCCGAGGCGATGGATTTGATTC
TCGTGGCAAGCGTGAATTTGATAGGCATAGTGAAGTGATAGATCTGGCCTGAAGCATGAGGACAAACGC
GGAGGTAGCGGCTCTCACACTGGGGAAGTGTCAAAGATGAATTAAGTATTTGGATCAATCAATGTGA
CTGAGGAAACACCTGAAGGTGAAGAGCACCTGTGGCAGATACTGAAAATAAGGAGAACGAAGTTGAAGA
GGTTAAGGAAGAGGGTCCAAAAGAGATGACTCTGGATGAGTGGAAAGCTATTCAAATAAAGACCGAGCA
AAAGTAGAATTTAATATCCGAAAACCAATGAAGCGCCGATGGACAATGAAAAAGGGATTTGTCTGC
ATAAATCAAAAAGTGAAGAGGCTCATGTGAAGATTGAGTATGAGCATCATTTCCGGAAGCCAGCAAA
TGATATAACATCTCAACTGGAGATCAATTTTGGAGACTTAGGCCGCCAGGACGTGGTGGCAGAGGAGGA
CGTGGTGGCGTGGGCGTGGTGGACGTCCAAACCGTGGCAGCAGGACTGATAAGTCAAGTGCTTCTGCTC
CTGATGTAGATGACCCAGAGGCTTCCAGCTCTGGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



View online »

Protein Sequence:	>MR206038 representing BC013665 Red=Cloning site Green=Tags(s) MPGHLQEGFGCVVTNRFQDLFDDSDPFEVLKAAENKKKEAGGGVGPGAKSAAQAAQTSNAAGKQL RKEQKDRKNPLPPSVGVADKKEETQPPVALKKEGIRRVGRRPDQQLQGDGKLIDRRAERRPPRERRFEK PLEEKEGEGEFSVDRPIIERPIRGRGGLGRGRGGRGMGRGDGFSRGKREFDRHSGSDRSGLKHEDKR GGSGSHNWGTVKDELTDLDQSNVTEETPEGEEHPVADTENKENEVEEVKEEGPKEMTLDEWKAIQNKDRA KVEFNIRKPNEGADGQWKGFVLHKSKEEAHAEDSVMDDHFRKPANDITSQLEINFGDLGRPGRGGRGG RGGRGRGGRPNRGSRTDKSSASAPDVEDDPEAFPALA TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	BC013665
ORF Size:	1158 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	BC013665.1
RefSeq Size:	1338 bp
RefSeq ORF:	1160 bp
Locus ID:	66870
Cytogenetics:	6 C1
MW:	49.1 kDa