

Product datasheet for **MR206068A1V**

Mouse Cct6a (BC018459) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Cct6a (BC018459) AAV Particle
Tag:	Myc-DDK
Symbol:	Cct6a
Synonyms:	CCT-zeta; CCT-zeta-1; Cct6; Cctz-1; TCP-1-zeta
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR206068 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGACAGAGAAACACTCATCGATGTGGCCAGGACATCTCTGCGGACTAAAGTTCATGCTGAACCTGCAG
ATGTCTTGACAGAGGCTGTAGTGGACTCCATCTTGCCATTAGGAAAAAGGACGAGCCATTGACCTCTT
CATGGTTGAGATCATGGAGATGAAGCATAAATCTGAGACAGATACAAGCTTAATCAGAGGGCTTGTTTTG
GATCATGGAGCTCGGCATCCTGATATGAAGAAGAGAGTGGAAAATGCCTACATCCTCACGTGCAACGTGT
CCTTAGAGTATGAGAAAACAGAAGTGAATTCTGGGTTTTTTTACAAGAGTGCAGAAGAGAGAGAAAAACT
AGTAAAGGCTGAAAGAAAATTCATTGAAGATAGAGTTAAAAAATCATAGAGCTGAAAAAGAAAGTCTGT
GGTGAATCAGATAAAGGATTTGTCGTTATTAATCAAAAGGGGATTGACCCCTTTTCCTTAGATGCCCTTG
CGAAAGAAGGGATCGTAGCTCTGCGCAGAGCCAAGAGGAGAAACATGGAGAGGCTGACACTTGCTTGTTGG
TGGGATAGCTCTGAATTCCTTTGATGACCTGAATCCTGACTGTTTGGGACATGCAGGGCTTGCTATGAG
TATACACTGGGTGAGGAGAAGTTCACCTTTATTGAGAAGTGAACAATCCCCGTTCTGTCACTTTACTGG
TTAAAGGACCAAATAAGCACACACTGACTCAAATCAAGGATGCAATAAGAGATGGCTTGAGGGCTGTCAA
AAATGCTATTGATGATGGCTGTGTTGTCCAGGTGCGGGTGCAGTAGAAGTGGCACTGGCAGAAGCTCTG
ATTAATACAAGCCAGTGTGAAGGGCAGGGCGCAGCTTGGAGTCCAGGCATTTGCAGATGCCTTGCTCA
TCATTCCAAGGTTCTTGCGCAAACTCTGGTTTTGACCTTCAGGAAACATTAGTTAAAGTTCAAGCTGA
ACATTCAGAATCGGGCCAGCTCGTAGGTGTGGATCTGAGCACAGGTGAGCCGATGGTGGCCGAGAGATG
GGTGTGTGGGATAACTACTGTGTGAAGAAGCAGCTGCTACACTCCTGTACTGTGATCGCCACCAACATTC
TCCTGGTCGACGAGATCATGCGAGCTGGGATGTCCTCTCTGAAGGGT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR206068 protein sequence Red=Cloning site Green=Tags(s) MDRETLIDVARTSLRTKVHAELADVL TEAVVDSILAIRKKDEPIDLFMVEIMEMKHKSETDTSIRGLVL DHGARHPDMKKRVENAYILT CNVSLEYEKTEVNSGFFYKSAEEREKL VKAERKFIEDRVKKI IELKKKVC GDSKGFVVINQKGIDPFSLDALAKEGIVALRRAKRRNMERL TLACGGIALNSFDDLNPDCLGHAGLVYE YTLGEEKFTFIEKCNNPRSVTLLVKGPKNKHTLTQIKDAIRDGLRAVKNAIDGCVVPGAGAVEVALAEAL IKYKPSVKGRAQLGVQAFADALLIIPKVL AQNSGFDLQETLVKVQAEHSESGQLVGVDLSTGEPMVAEM GVWDNYCVKKQLLHSCCTVIATNILLVDEIMRAGMSSLKG TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	BC018459
ORF Size:	1167 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>BC018459</u> , <u>AAH18459</u>
RefSeq Size:	1795 bp
RefSeq ORF:	1169 bp
Locus ID:	12466
Cytogenetics:	5 G1.3
MW:	42.9 kDa