

Product datasheet for **MR206068L3V**

Cct6a (BC018459) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cct6a (BC018459) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cct6a
Synonyms:	CCT-zeta; CCT-zeta-1; Cct6; Cctz-1; TCP-1-zeta
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC018459
ORF Size:	1167 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206068).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC018459 , AAH18459
RefSeq Size:	1795 bp
RefSeq ORF:	1169 bp
Locus ID:	12466
Cytogenetics:	5 G1.3
Gene Summary:	Component of the chaperonin-containing T-complex (TRiC), a molecular chaperone complex that assists the folding of proteins upon ATP hydrolysis. The TRiC complex mediates the folding of WRAP53/TCAB1, thereby regulating telomere maintenance. The TRiC complex plays a role in the folding of actin and tubulin.[UniProtKB/Swiss-Prot Function]


[View online »](#)