

Product datasheet for **MR206081L4V**

Pafah2 (NM_133880) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pafah2 (NM_133880) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pafah2
Synonyms:	2310074E22Rik; AI507170
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_133880
ORF Size:	1173 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206081).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_133880.1 , NP_598641.1
RefSeq Size:	3183 bp
RefSeq ORF:	1173 bp
Locus ID:	100163
UniProt ID:	Q8VDG7
Cytogenetics:	4 D3
Gene Summary:	Has a marked selectivity for phospholipids with short acyl chains at the sn-2 position. May share a common physiologic function with the plasma-type enzyme (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)