

Product datasheet for **MR206095L3V**

Usp50 (NM_029163) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Usp50 (NM_029163) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Usp50
Synonyms:	1700086G18Rik; 4930511O11Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_029163
ORF Size:	1173 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206095).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_029163.1 , NP_083439.1 , NP_083439.2
RefSeq Size:	1748 bp
RefSeq ORF:	1173 bp
Locus ID:	75083
UniProt ID:	Q6P8X6
Cytogenetics:	2 F1
Gene Summary:	May recognize and hydrolyze the peptide bond at the C-terminal Gly of ubiquitin. Involved in the processing of poly-ubiquitin precursors as well as that of ubiquitinated proteins (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)