

Product datasheet for **MR206101L4V**

Dpf2 (NM_011262) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dpf2 (NM_011262) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dpf2
Synonyms:	2210010M07Rik; BAF45D; Req; ubi-d4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_011262
ORF Size:	1176 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206101).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011262.2
RefSeq Size:	3314 bp
RefSeq ORF:	1176 bp
Locus ID:	19708
UniProt ID:	Q61103
Cytogenetics:	19 4.34 cM



[View online »](#)

Gene Summary:

Plays an active role in transcriptional regulation by binding modified histones H3 and H4. Is a negative regulator of myeloid differentiation of hematopoietic progenitor cells (By similarity). Might also have a role in the development and maturation of lymphoid cells (PubMed:7961935). Involved in the regulation of non-canonical NF-kappa-B pathway (By similarity).[UniProtKB/Swiss-Prot Function]