

Product datasheet for **MR206131L3V**

Golm1 (NM_027307) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Golm1 (NM_027307) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Golm1
Synonyms:	2310001L02Rik; AW125446; D030064E01Rik; Golph2; GP73; PSEC0257
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_027307
ORF Size:	1182 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR206131).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_027307.4
RefSeq Size:	3693 bp
RefSeq ORF:	1182 bp
Locus ID:	105348
UniProt ID:	Q91XA2
Cytogenetics:	13 B2
Gene Summary:	Unknown. Cellular response protein to viral infection (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)